



Certificate

The following products:

DUPLEX IV / DUPLEX IV PROFESSIONAL
ELO-RAYEX / HF-RAYEX / CAR-RAYEX
MINI-RAYONEX / MINI-RAYONEX 5G
E-SMOG-RAYONATOR II / HF-RAYONATOR III / 5G-RAYONATOR / CHAKRA-RAYONATOR

manufactured by

RAYONEX BIOMEDICAL GmbH

Sauerland-Pyramiden 1
57368 Lennestadt
GERMANY

were accurately examined in our institute regarding safety and efficacy.

We hereby confirm that the products listed above meet the requirements for an energetic harmonization in the living organism according to the principles of Bioresonance according to Paul Schmidt.

The results of our examination are listed below by product:

DUPLEX IV / DUPLEX IV PROFESSIONAL

The frequency spectra used in the DUPLEX IV devices can reduce strains caused by geopathy, underground watercourses, dislocations, Curry-nets, Hartman-grids, and Benker-grids in the organism. We were able to confirm an efficiency range of up to 20 metres in diameter in the corresponding test setup.

ELO-RAYEX

In the ELO-RAYEX the six (6) most important fundamental frequency values are installed, able to harmonize the impact from low-frequency electromagnetic pollution in the organism. We were able to confirm an efficiency range of up to 2 meters in diameter in the corresponding test setup.

HF-RAYEX

In the HF-RAYEX the six (6) most important fundamental frequency values are installed, able to harmonize the impact caused particularly from high-frequency pulse modulated influences in the organism. We were able to confirm an efficiency range of up to 2 meters in diameter in the corresponding test setup.

CAR-RAYEX

In the CAR-RAYEX the six (6) most important fundamental frequency values are installed, able to harmonize the impact caused particularly from Electro-smog (3) and pulse modulated (3) influences in the organism. We were able to confirm an efficiency range of up to 2 meters in diameter in the corresponding test setup.

MINI-RAYONEX

In the MINI-RAYONEX the most important fundamental frequency value (12,5) for the stimulation of the energy household is installed. As each strain of the organism needs energy, the MINI-RAYONEX can be used in a very various way. We were able to confirm an efficiency range of up to 2,5 meters in diameter in the corresponding test setup.

MINI-RAYONEX 5G

The MINI-RAYONEX 5G builds on the previous MINI-RAYONEX and its fundamental frequency value of 12,5. In addition, however, it contains a further fundamental frequency value, which in the experience is always required when you are under the influence of a 5G network.

The MINI-RAYONEX 5G therefore contains two fundamental frequency values of Bioresonance according to Paul Schmidt: firstly, the fundamental frequency value 12,5, and secondly the most important 5G fundamental frequency value of Bioresonance according to Paul Schmidt. We were able to confirm an efficiency range of up to 2,5 meters in diameter in the corresponding test setup.

E-SMOG-RAYONATOR II

The E-SMOG-RAYONATOR II contains, as the ELO-RAYEX does, the six (6) most important fundamental frequency values referring to the impact of low-frequency electromagnetic pollution. Contrary to the ELO-RAYEX, the E-SMOG-RAYONATOR II delivers its frequencies via a connection to the DUPLEX IV / DUPLEX IV PROFESSIONAL. In the corresponding test set-up, we were able to confirm an efficiency range equal to DUPLEX IV / DUPLEX IV PROFESSIONAL.

HF-RAYONATOR III

The HF-RAYONATOR III contains 10 important fundamental frequency values related to the pulse modulated impact (like Wi-Fi, Bluetooth, Mobile phones etc.). Contrary to the HF-RAYEX, the HF-RAYONATOR III delivers its frequencies via a connection to the DUPLEX IV / DUPLEX IV PROFESSIONAL. In the corresponding test set-up, we were able to confirm an efficiency range equal to DUPLEX IV / DUPLEX IV PROFESSIONAL.

5G-RAYONATOR

The 5G-RAYONATOR contains 10 dipole antenna systems of Bioresonance according to Paul Schmidt inside the device. Eight of them are related to the specific impact of 5G radiation. The other two are the fundamental frequency values 12,50 and 21,50. The 5G-RAYONATOR delivers its frequencies via a connection to the DUPLEX IV / DUPLEX IV PROFESSIONAL. In the corresponding test set-up, we were able to confirm an efficiency range equal to DUPLEX IV / DUPLEX IV PROFESSIONAL.

CHAKRA-RAYONATOR

The CHAKRA-RAYONATOR contains a total of 10 dipole antenna systems of Bioresonance according to Paul Schmidt inside the device. Seven of them each refer to the fundamental frequency value of one of the chakras. The other three contain the fundamental frequency values 10,00; 12,50 and 21,50. The CHAKRA-RAYONATOR delivers its frequencies via a connection to the DUPLEX IV / DUPLEX IV PROFESSIONAL. In the corresponding test set-up, we were able to confirm an efficiency range equal to DUPLEX IV / DUPLEX IV PROFESSIONAL.



Head of Institute
Institut für Bioenergetik | Bioenergetik SH Nübbel

Wolfgang Sievers, 17.11.2020